Medical microbiology

Lab 7

Bacterial motility



Dr. Sarwa Aziz Khalid

Assis. Lec . Shahad Muaad Tawfeeq

Motility

Motility occurs in rod shaped bacteria and spirochetes but is almost never found in cocci . The major organelles of motility in bacteria are flagella . Flagella allow cells to move toward nutrients in the environment or move away from harmful substances , such as acids , in a complicated process called chemotaxis .

Motility can be determined by several methods:

1-Wet Mount method

2-Hanging Drop

technique 3-Semisolid


media

Most bacterial microscopic preparations result in death of the microorganisms due to heat-fixing and staining. Simple wet mounts and the hanging drop technique allow observation of living cells to determine motility.

A wet mount preparation is made by placing the specimen in a drop of water on a microscope slide and covering it with a cover glass then examin it under microscope.


A hanging drop preparation allows longer observation of the specimen since it doesn't dry out as quickly. A thin ring of petroleum jelly is applied around the well of a depression slide. A drop of water is then placed in the center of the cover glass and living microbes are transferred into it. The depression microscope slide is carefully placed over the cover glass in such a way that the drop is received into the depression and is undisturbed. The petroleum jelly causes the cover glass to stick to the slide.

The preparation may then be picked up, inverted so the cover glass is on top, and placed under the microscope for examination



Hanging drop method

Semisolid media: Inoculate tubes of semisolid media (contain agar about 0.3-0.5 %) with the organisms by stab (using needle) and incubate it for 24 hours.



Left positive and right negative